

## Instructions for Use



### STEP 1

Find the lower section of the IPC (marked by the fluorescent indicator). While ensuring the large buckle remains to the right of the victim, raise the victim's knees and place the IPC perpendicular to the body. The fluorescent indicator should be in the center next to the groin area.



### STEP 2

Close both leg straps around the victim's legs and slide the upper part of the IPC carefully under victim's back.



### STEP 3

Reach into the under arm areas and pull each side of the main harness around the victim's body, creating shoulder straps.



### STEP 4

Take a seat between the victim's legs and place arms into each of the shoulder straps. Get your shoulders as close to the victim's shoulders as possible.

Make sure the shoulder straps are snug (like backpack straps), tightening the large buckle if necessary.

For optimal comfort, make sure the shoulder pads are correctly located and close the chest buckle.



### STEP 5

Pulling the victim's arms close to you, roll over onto all four (the victim will be on your back). Slowly transfer weight to your legs and stand up, keep your back as straight as possible.



# INJURED PERSONNEL CARRIER IPC™



The IPC allows you to carry an injured fellow hiker long distances in a way that is more comfortable than any existing method, both for the carrier and the evacuee.

- The only evacuation method that keeps both hands free
- Strong yet extremely compact and lightweight
- Ideal for soldiers, hikers, fire fighters and emergency responder

## MEDIHELP HEALTHCARE PVT. LTD.

D-6, 6052, Vasant Kunj, New Delhi-110070

Ph.: +91-11-4605 2020 Mobile: +91 9871 555 932

E-mail: medihelpindia@gmail.com Website: www.medihelpindia.in

