

I-gel

A-06.



 **i-gel**
PLUS
i-gel® has evolved





15mm connector

Reliable connection to any standard catheter mount or connection

Colour coded hook ring

To secure the i-gel® Plus in position with the airway support strap

Proximal end of gastric channel

Clearly displayed product information

For quick easy reference. Includes confirmation of size and weight guidance

Position guide

Supplementary oxygen port

For the administration of passive oxygenation as a component of cardiocerebral resuscitation (CCR)

Gastric channel

Enhances patient safety by providing a mechanism for the management of regurgitant fluid

Integral bite block

Reduces the possibility of airway channel occlusion

Buccal cavity stabiliser

Aids insertion and eliminates the potential for rotation

Larger gastric channel

To enable insertion of a larger gastric tube

+ i-gel® Plus added benefit

Intubation ramp

To optimise use as a conduit for intubation

+ i-gel® Plus added benefit

Longer cuff tip

To help enhance the oesophageal seal

+ i-gel® Plus added benefit

Epiglottic rest

Reduces the possibility of epiglottis 'down folding' and airway obstruction

Non-inflatable cuff

Eliminates the need for cuff inflation after insertion, allowing easy and rapid insertion

Distal end of gastric channel

